

# SAPPHIRE<sup>®</sup> II PRO

Coronary Dilatation Catheter

Cross with Confidence  
Every Time

1.0mm-1.5mm



Robust  
Material

Superb  
Balloon  
Memory


Ultimate  
Crossability

*Redefining the Standard of Balloon Angioplasty*


# SAPPHIRE<sup>®</sup> II PRO

Coronary Dilatation Catheter

## Tip Entry Comparison\*



## Crossing Profile Comparison\*



\*Data on file, scale depictions

## Ordering Information


Ø (mm)	Length (mm)					
	5	8	10	12	15	20
1.0	210-053-5UJ	210-083-5UJ	210-103-5UJ	—	210-153-5UJ	—
1.25	212-053-5UJ	212-083-5UJ	212-103-5UJ	—	212-153-5UJ	—
1.5	—	—	215-103-5UJ	215-123-5UJ	215-153-5UJ	215-203-5UJ

## Compliance Chart

Pressure (atm)	Ø (mm)		
	1.0	1.25	1.5
2	0.94	1.18	1.42
4	0.97	1.22	1.46
6 NOM*	1.00	1.25	1.50
8	1.03	1.28	1.54
10	1.06	1.32	1.58
12	1.09	1.35	1.62
14	1.12	1.38	1.66
16 RBP**	1.15	1.42	1.70
18	1.18	1.45	1.74
20	1.21	1.48	1.78
22	1.24	1.52	1.82

\* Nominal Pressure. The nominal in-vitro device specifications do not take into account any lesion resistance.

\*\* Rated Burst Pressure. Do not exceed RBP.



For more information please visit our website at [www.OrbusNeich.com](http://www.OrbusNeich.com) or contact us:

### CORPORATE HEADQUARTERS

13th Floor, Harcourt House  
39 Gloucester Road  
Wanchai, Hong Kong  
Phone +852.2802.2288  
Fax +852.2507.3532

### EMEA HEADQUARTERS

Drs. W. van Royenstraat 5  
3871 AN Hoevelaken  
The Netherlands  
Phone +31.33.254.1150  
Fax +31.33.254.1151  
Toll Free Phone 00800.0254.1150\*  
Toll Free Fax 00800.0254.1151\*

\* Belgium, Denmark, France, Germany, Ireland, Netherlands, Norway, Sweden and UK

© 2018 OrbusNeich. All rights reserved.

#G-70-0410 Rev03

Sapphire II PRO and OrbusNeich are Trademarks of OrbusNeich Medical, Inc.

 **OrbusNeich<sup>®</sup>**  
Pioneers in life-changing technologies