

Ordering Information




Balloon Diameter (mm)		Balloon Working Length (mm)						
		5	8	10	12	15	20	30
OTW	1.0	501000512	501000812	501001012	-	501001512	-	-
Rx	1.0	210-053-5UU	210-083-5UU	210-103-5UU	-	-	-	-
	1.25	212-053-5UU	212-083-5UU	212-103-5UU	-	212-153-5UU	-	-
	1.5	-	-	215-103-5UU	215-123-5UU	215-153-5UU	215-203-5UU	-
	1.75	-	-	217-103-5UU	-	217-153-5UU	217-203-5UU	-
	2.0	-	-	220-103-5UU	220-123-5UU	220-153-5UU	220-203-5UU	-
	2.25	-	-	222-103-5UU	-	222-153-5UU	222-203-5UU	-
	2.5	-	-	225-103-5UU	225-123-5UU	225-153-5UU	225-203-5UU	225-303-5UU
	2.75	-	-	227-103-5UU	-	227-153-5UU	227-203-5UU	-
	3.0	-	-	230-103-5UU	230-123-5UU	230-153-5UU	230-203-5UU	230-303-5UU
	3.25	-	-	232-103-5UU	-	232-153-5UU	232-203-5UU	-
	3.5	-	-	235-103-5UU	-	235-153-5UU	235-203-5UU	235-303-5UU
	4.0	-	-	240-103-5UU	-	240-153-5UU	240-203-5UU	-

Compliance Chart


Pressure (atm)	Balloon Diameter (mm)											
	1.0	1.25	1.5	1.75	2.0	2.25	2.5	2.75	3.0	3.25	3.5	4.0
2	0.83	1.04	1.29	1.41	1.66	1.89	2.12	2.33	2.50	2.74	2.98	3.38
4	0.94	1.17	1.42	1.62	1.87	2.12	2.36	2.60	2.82	3.06	3.31	3.77
6 NOM*	1.00	1.25	1.50	1.75	2.00	2.25	2.50	2.75	3.00	3.25	3.50	4.00
8	1.04	1.30	1.56	1.84	2.09	2.34	2.60	2.86	3.13	3.38	3.64	4.16
10	1.08	1.35	1.60	1.91	2.16	2.42	2.68	2.94	3.23	3.49	3.74	4.29
12	1.11	1.38	1.63	1.96	2.21	2.48	2.74	3.01	3.31	3.57	3.83	4.39
14 RBP**	1.13	1.41	1.66	2.01	2.26	2.53	2.79	3.07	3.38	3.64	3.90	4.48
16	1.15	1.41	1.69	2.05	2.30	2.57	2.84	3.12	3.44	3.70	3.96	4.55
18	1.17	1.46	1.71	2.09	2.34	2.61	2.88	3.17	3.50	3.76	4.02	4.62

* Nominal Pressure. The nominal in-vitro device specifications do not take into account any lesion resistance.
** Rated Burst Pressure. Do not exceed RBP.

Over the Wire (OTW)



Rapid Exchange (Rx)



For more information, please contact your local CSI representative or call customer service at 1-877-274-0901.

Caution: Federal Law (USA) restricts the sale of these devices to sale by, or on the order of, a physician.

CSI is a registered trademark of Cardiovascular Systems, Inc. ©2020 Cardiovascular Systems, Inc. All rights reserved. CSI EN-4771.G 0620 Sapphire II PRO and OrbusNeich are registered trademarks of OrbusNeich Medical Corporation. ©2020 OrbusNeich. G-70-1125. All other trademarks cited herein are the trademarks of their respective owners.

Distributed by:

CSI
CARDIOVASCULAR
SYSTEMS, INC.

1225 Old Hwy 8 NW
St. Paul, MN 55112
T: 651.259.1600
F: 612.677.3355
www.csi360.com

Manufactured by:

OrbusNeich
Pioneers in life-changing technologies

1 Jinkui Road,
Fuitian Free Trade Zone,
Shenzhen, 518038, China
T: +86.755.8358.0181
F: +86.755.8358.0169
www.OrbusNeich.com



Redefining the Standard of
Balloon Angioplasty

Distributed by:

CSI
CARDIOVASCULAR
SYSTEMS, INC.

1 Jinkui Road,
Fuitian Free Trade Zone,
Shenzhen, 518038, China
T: +86.755.8358.0181
F: +86.755.8358.0169
www.OrbusNeich.com

Manufactured by:

OrbusNeich
Pioneers in life-changing technologies


1.0 mm - 4.0 mm Semi-compliant Balloon for Vascular Treatment

Cross with Confidence

Specifically engineered for crossing difficult lesions and tracking tortuous anatomy.


P-TECH Shaft for Optimal Pushability and Kink Resistance

The P-TECH Shaft is uniquely constructed with a small catheter profile that is kink resistant and facilitates excellent deliverability.



Sub-zero Tip Technology for Smooth Lesion Entry

The well balanced sub-zero tapered tip has an ultra low profile and provides smooth entry through the tightest lesions.




Smallest Crossing Profile in the US Market

XR Balloon Engineered for Crossability

The proprietary "XR balloon" facilitates best-in-class crossability and recrossability without compromising durability and robustness.

Tip Entry and Crossing Profile Comparison



OrbusNeich Data on file.