



Ordering Information

014

Balloon Diameter (mm)	Catheter Working Length – 150 cm								Sheath Compatibility
	20	40	60	80	100	120	150	240	
1.5	581502022	581504022	-	581508022	-	581512022	-	581524022	4F
2.0	582002022	582004022	582006022	582008022	582010022	582012022	582015022	582024022	
2.5	582502022	582504022	582506022	582508022	582510022	582512022	582515022	582524022	
3.0	583002022	583004022	583006022	583008022	583010022	583012022	583015022	583024022	
3.5	583502022	583504022	583506022	583508022	583510022	583512022	583515022	583524022	
4.0	584002022	584004022	584006022	584008022	584010022	584012022	584015022	584024022	
5.0	585002022	585004022	585006022	585008022	585010022	585012022	585015022	585024022	5F
6.0	586002022	586004022	586006022	586008022	586010022	586012022	586015022	586024022	6F
Balloon Diameter (mm)	Catheter Working Length – 200 cm								Sheath Compatibility
	20	40	60	80	100	120	150	240	
2.0	582002032	582004032	582006032	582008032	582010032	582012032	582015032	582024032	4F
2.5	582502032	582504032	582506032	582508032	582510032	582512032	582515032	582524032	
3.0	583002032	583004032	583006032	583008032	583010032	583012032	583015032	583024032	
3.5	583502032	583504032	583506032	583508032	583510032	583512032	583515032	583524032	
4.0	584002032	584004032	584006032	584008032	584010032	584012032	584015032	584024032	
5.0	585002032	585004032	585006032	585008032	585010032	585012032	585015032	585024032	5F
6.0	586002032	586004032	586006032	586008032	586010032	586012032	586015032	586024032	6F

018

Balloon Diameter (mm)	Catheter Working Length – 150 cm								Sheath Compatibility
	20	40	60	80	100	120	150	240	
2.0	592002022	592004022	592006022	592008022	592010022	592012022	592015022	592024022	4F
2.5	592502022	592504022	592506022	592508022	592510022	592512022	592515022	592524022	
3.0	593002022	593004022	593006022	593008022	593010022	593012022	593015022	593024022	
3.5	593502022	593504022	593506022	593508022	593510022	593512022	593515022	593524022	
4.0	594002022	594004022	594006022	594008022	594010022	594012022	594015022	594024022	
5.0	595002022	595004022	595006022	595008022	595010022	595012022	595015022	595024022	5F
6.0	596002022	596004022	596006022	596008022	596010022	596012022	596015022	596024022	
7.0	597002022	597004022	597006022	597008022	597010022	-	597015022	-	6F
Balloon Diameter (mm)	Catheter Working Length – 200 cm								Sheath Compatibility
	20	40	60	80	100	120	150	240	
2.0	592002032	592004032	592006032	592008032	592010032	592012032	592015032	592024032	4F
2.5	592502032	592504032	592506032	592508032	592510032	592512032	592515032	592524032	
3.0	593002032	593004032	593006032	593008032	593010032	593012032	593015032	593024032	
3.5	593502032	593504032	593506032	593508032	593510032	593512032	593515032	593524032	
4.0	594002032	594004032	594006032	594008032	594010032	594012032	594015032	594024032	
5.0	595002032	595004032	595006032	595008032	595010032	595012032	595015032	595024032	5F
6.0	596002032	596004032	596006032	596008032	596010032	596012032	596015032	596024032	
7.0	597002032	597004032	597006032	597008032	597010032	-	597015032	-	6F

035

Balloon Diameter (mm)	Catheter Working Length – 135 cm							Sheath Compatibility
	20	40	60	80	100	150	240	
3.0	603002022	603004022	603006022	603008022	-	603015022	-	5F
4.0	604002022	604004022	604006022	604008022	604010022	604015022	604024022	
5.0	605002022	605004022	605006022	605008022	605010022	605015022	605024022	6F
6.0	606002022	606004022	606006022	606008022	606010022	606015022	606024022	
7.0	607002022	607004022	607006022	607008022	607010022	607015022	-	7F
Balloon Diameter (mm)	Catheter Working Length – 200 cm							Sheath Compatibility
	20	40	60	80	100	150	240	
5.0	605002032	605004032	605006032	605008032	605010032	605015032	605024032	6F
6.0	606002032	606004032	606006032	606008032	606010032	606015032	606024032	
7.0	607002032	607004032	607006032	607008032	607010032	607015032	-	7F

Strength and non-compliant precision
without compromising deliverability

For more information, please contact your local CSI representative or call customer service at 1-877-274-0901.

Caution: Federal Law (USA) restricts the sale of these devices to sale by, or on the order of, a physician.

CSI is a registered trademark of Cardiovascular Systems, Inc. ©2021 Cardiovascular Systems, Inc. All rights reserved. CSI EN-6547.A 0421. JADE and OrbusNeich are registered trademarks of OrbusNeich Medical Company Limited. ©2021 OrbusNeich Medical Company Limited or its affiliates. G-70-1334 Rev02. All other trademarks cited herein are the trademarks of their respective owners.

Distributed by:



CARDIOVASCULAR
SYSTEMS, INC.

1225 Old Hwy 8 NW
St. Paul, MN 55112

T: 651.259.1600
877.274.0901

F: 612.677.3355
www.csi360.com

Manufactured by:



OrbusNeich®
Pioneers in life-changing technologies

1 Jinkui Road,
Futian Free Trade Zone,
Shenzhen, 518038, China

T: +86.755.8358.0181
F: +86.755.8358.0169
www.OrbusNeich.com

Distributed by:




CARDIOVASCULAR
SYSTEMS, INC.

CARDIOVASCULAR
SYSTEMS, INC.

Manufactured by:




OrbusNeich®
Pioneers in life-changing technologies



The **JADE PTA** balloon catheter series is an over-the-wire platform built on strength and non-compliant precision without compromising deliverability during peripheral interventions.

JADE PTA balloon series offers a **Diverse Balloon Portfolio**


JADE is the only non-compliant balloon platform available in all three guidewire compatible systems on the US market



Available in **200 cm** shaft lengths for radial access and balloon working lengths from **20 - 240 mm**



Designed for **Excellent Deliverability**

Z-Tip technology offers exceptional entry and competitive crossing profile while proprietary dual balloon coating provides enhanced trackability



Non-Compliant Balloon for **Safety and Accuracy**

Consistent uniform pressure distribution (minimizing dog-boning effect) during balloon dilatation in all plaque morphologies



Up to **20 atm RBP** for opening resistant lesions with quick deflation time

Z-Tip technology for crossing challenging lesions

Compliance Chart

Pressure (atm)	Balloon Diameter (mm)								
	1.5 ¹	2.0	2.5	3.0	3.5 ²	4.0	5.0	6.0	7.0 ³
2	1.41	1.88	2.37	2.82	3.27	3.70	4.67	5.60	6.64
4	1.43	1.90	2.40	2.86	3.32	3.76	4.74	5.68	6.73
6	1.45	1.93	2.42	2.89	3.36	3.82	4.80	5.76	6.82
8	1.47	1.95	2.45	2.93	3.41	3.88	4.87	5.84	6.91
10 NOM*	1.48	1.98	2.47	2.96	3.45	3.94	4.93	5.92	7.00
12 NOM*	1.50	2.00	2.50	3.00	3.50	4.00	5.00	6.00	7.09
14	1.52	2.02	2.53	3.04	3.55	4.06	5.07	6.08	7.18
16 RBP**	1.53	2.05	2.55	3.07	3.59	4.12	5.13	6.16	7.27
18 RBP**	1.55	2.07	2.58	3.11	3.64	4.18	5.20	6.24	7.36
20 RBP**	1.57	2.10	2.60	3.14	3.68	4.24	5.26	6.32	7.45
22	1.59	2.12	2.63	3.18	3.73	4.30	5.33	6.40	-
24	1.60	2.14	2.65	3.22	3.78	4.36	5.40	-	-
26	1.62	2.17	2.68	3.25	3.82	4.42	-	-	-

¹ 1.5 mm configurations available only for 014 platform
² 3.5 mm configurations available only for 014 and 018 platforms
³ 7.0 mm configurations available only for 018 and 035 platforms
* Nominal Pressure. The nominal in-vitro device specifications do not take into account any lesion resistance.
** Rated Burst Pressure. Do not exceed RBP.